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Short instruction

1. Start of an existing program for the TwinCubator® with amplification block

- 1) switch on the instrument
- 2) wait for initialization of the instrument (see paragraph 3)
- 3) in order to access a program select the respective program with help of the **+** and **-** keys and press the Start/Stop key
- 4) the instrument switches automatically from one program step to the next
- 5) an acoustic signal will indicate the end of the program; the instrument will stop shaking

2. Start of an existing program for the TwinCubator® with hybridization block

- 1) switch on the instrument
- 2) wait for initialization of the instrument (see paragraph 4)
- 3) in order to access a program select the respective program with help of the **+** and **-** keys and press the Start/Stop key
- 4) an acoustic signal indicates the end of an incubation step; shaking and heating will be continued
- 5) by pressing the **→** key the shaking function will be stopped and the program will proceed to the next program step
- 6) by pressing the **→** key again the next program step will be started
- 7) repeat step 4 to 6 until the last program step is reached
- 8) an acoustic signal will indicate the end of the program; the instrument will stop shaking

3. Initialization of TwinCubator® with amplification block comprises the 4 steps:

- 1) the red and yellow LED lights of the Start/Stop key and the programming key, respectively, will come on
- 2) an acoustic signal sounds
- 3) display indicates "Twin Sarstedt Block Vx.yy WW/YY"
- 4) display switches over to basic setting (upper line: indication of set and actual temperature, time and shaking frequency; lower line: indication of name of step, step mode, program number and program step number)

4. Initialization of TwinCubator® with hybridization block comprises the 4 steps:

- 1) the red and yellow LED lights of the Start/Stop key and the programming key, respectively, will come on
- 2) an acoustic signal sounds
- 3) display indicates "TwinCubator Waennchen Vx.yy WW/JJ"
- 4) display switches over to basic setting (upper line: indication of set and actual temperature, time and shaking frequency; lower line: indication of name of step, step mode, program number and program step number)