GenoType® DNA Isolation Kit

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GenoType® DNA Isolation Kit **Kit for isolation and purification of genomic DNA from**

- human whole blood**
- body fluids**
- tissue**
- eukaryotic cell cultures**

Methodology

The **GenoType® DNA Isolation Kit** permits the easy and fast manual isolation and purification of genomic DNA from human whole blood, body fluids, tissue, and eukaryotic cell culture material for *in vitro* diagnostic analysis in 10-20 min (after lysis). Human genomic DNA can be purified from fresh or frozen EDTA or citrate (not heparin) whole blood from common blood collection systems. Additionally, frozen tissue samples, cell suspensions, cell pellets, and body fluids like amniotic fluid, synovial fluid etc. may be used.

This kit must not be used for DNA isolation from stool samples. It was neither validated for isolation from bacteria, fungi, or viruses nor for isolation of RNA.

The whole procedure is divided into four steps: 1. lysis of cells under anti-chaotropic conditions at elevated temperature, 2. binding of genomic DNA to the membrane of the Spin Filter, 3. washing of the membrane to remove contaminations and ethanol removal, 4. elution of the purified genomic DNA. Due to its purification grade, the isolated genomic DNA may directly be used for a multitude of downstream applications or can be stored at –20°C.

Storage and Precautions

Store all kit components at room temperature (18-25°C). Close ethanol-containing Wash Buffer I and II appropriately and store at room temperature. Store reconstituted Proteinase K at –20°C. In order to avoid repeated freezing and thawing, aliquot Proteinase K.

Patient specimens must always be considered as potentially infectious. Samples from risk patients must always be labeled and handled under suitable safety conditions. Observe all federal, state, and local safety and environmental regulations. Always wear suitable protective clothing and gloves.

Observe the usual precautions for nucleic acid applications. It is essential that all reagents and materials used for DNA isolation are free from DNases.

When handling kit reagents, the following special safety measures must be applied:

Lysis Buffer M contains Ammonium Chloride and is harmful (R36 and S2-24)

Binding Buffer B6 contains 2-Propanol and is irritating and highly flammable (R11-36-67 and S2-7-16-S24/25/26).

Proteinase K is harmful (R36/37/38 and S2-22-24-36/37-42).

Wash Buffer I contains Guanidine Thiocyanat and is harmful (R20/21/22-32 and S2-13).

For additional information, please refer to material safety data sheets which can be downloaded from:

<http://www.hain-lifescience.de/products/msds.html>

Test Procedure

Preparation

Preheat thermal mixer to 56°C. Preheat the amount of Elution Buffer D needed (200 µl for each sample) to 56°C. Label the tubes needed: one 1.5 ml reaction tube (two for tissue protocol; not provided), one 1.5 ml Receiver Tube (white), and two 2.0 ml Receiver Tubes (yellow) for each sample processed. For each sample, place a Spin Filter in a labeled 2.0 ml Receiver Tube.

If the kit is used for the first time, reconstitute Proteinase K lyophilisate as follows: for the 12 tests kit, add 250 µl dd water to lyophilisate, for the 48 tests kit, add 550 µl dd water to lyophilisate, and vortex for approximately 5 sec. After use, store at –20°C. Thaw if already dissolved. Prepare wash buffers as follows: add 30 ml 99.8% ethanol p.a. to Wash Buffer I and 42 ml 99.8 ml ethanol p.a. to Wash Buffer II, respectively, and mix thoroughly. Store appropriately closed at room temperature.

The reagents must be free from precipitates; if necessary, carefully warm to up to max. 30°C to dissolve.

Notes

Discard pipette tips after each pipetting step.

To prevent contamination, change gloves if they have got in contact with sample material; open only one tube at a time; avoid wetting the rim of tubes with sample material.

Do not touch the Spin Filter membrane with pipette tip. During centrifugation steps, the Spin Filter is closed with the lid of the respective Receiver Tube.

All centrifugation steps are carried out at room temperature.

Using a lower volume of Elution Buffer D leads to a higher concentration of the purified DNA. Eluting twice each with 100 µl of Elution Buffer D may lead to a higher concentration of purified DNA.

For validating the correct use of this test kit a previously isolated DNA can be used as a positive control.

Prior to any subsequent application, the isolated DNA should be centrifuged for 1 min at full speed.

The reagents and Spin Filters of this kit are for single use only.

DNA Isolation from Whole Blood

Blood samples can be stored at room temperature for up to 6 h, or at 2-8°C for up to 24 h. For long time storage, freeze samples at –20°C or –80°C. Avoid repeated freezing and thawing. If cryoprecipitates appear after thawing, do not transfer them into the Spin Filter.

Various different primary tubes and anticoagulants (except for heparin) can be used to collect blood samples for use in this test kit. The amount of DNA that can be isolated with this kit from max. 50 µl of whole blood depends on the leukocyte count. The sample volume for the lysis step must be 50 µl; add water (molecular biology grade) if necessary.

Perform initial steps as indicated in section Preparation.

1. Transfer 1-50 μ l whole blood into a labeled 1.5 ml reaction tube (not provided); if necessary, add water (molecular biology grade) to a final volume of 50 μ l. Add 50 μ l Lysis Buffer M and 10 μ l Proteinase K (see Preparation) and vortex reaction tube thoroughly for at least 5 sec.
2. Incubate sample for 5 min at 56°C and 800 rpm in a thermal mixer.
3. Add 100 μ l Binding Buffer B6 and mix by pipetting up and down four times. Place a Spin Filter in a labeled 2.0 Receiver Tube (yellow) and pipette lysate on membrane of Spin Filter. Close lid and centrifuge Receiver Tube with Spin Filter in a standard table top centrifuge for 1 min at 13,000 x g.
4. Open lid and add 300 μ l Wash Buffer I (see Preparation) to membrane of Spin Filter. Close lid and centrifuge for 30 sec at 13,000 x g.
5. Open lid and discard Receiver Tube with filtrate. Place Spin Filter in a new, labeled 2.0 ml Receiver Tube and add 750 μ l Wash Buffer II (see Preparation). Close lid and centrifuge for 30 sec at 13,000 x g.
6. Open lid and discard filtrate. Place Spin Filter in the same Receiver Tube, close lid and centrifuge for 2 min at 13,000 x g to eliminate residual ethanol.
7. Transfer Spin Filter to a new, labeled 1.5 ml receiver tube (white). Add 200 μ l prewarmed Elution Buffer D to membrane of Spin Filter, close lid and incubate at room temperature for 1 min.
8. To elute DNA, centrifuge for 1 min at 6,000 x g.

DNA Isolation from Body Fluids or Eukaryotic Cell Suspensions

Best results with body fluids (e.g. synovial fluid, amniotic fluid) and cell suspensions are obtained from fresh cells. For long time storage, freeze material at -20°C or -80°C immediately after withdrawal or culture-harvesting, respectively. Avoid repeated freezing and thawing since this might accelerate the degradation process. Inappropriate storage of sample material leads to decreased DNA quantity and quality. The DNA amount that can be purified with this kit from 50 μ l body fluid or 10²-10⁵ cells depends on type and number of cells. The sample volume for the lysis step must be 50 μ l; add water (molecular biology grade) if necessary.

Perform initial steps as indicated in section Preparation.

1. Transfer max. 50 μ l body fluid or cell suspension into a labeled 1.5 ml reaction tube (not provided); if necessary, add water (molecular biology grade) to a final volume of 50 μ l. Add 50 μ l Lysis Buffer M and 10 μ l Proteinase K (see Preparation) and vortex reaction tube thoroughly for at least 5 sec.
2. Incubate sample for 5 min at 56°C and 800 rpm in a thermal mixer.
3. Add 100 μ l Binding Buffer B6 and mix by pipetting up and down four times. Place a Spin Filter in a labeled 2.0 Receiver Tube (yellow) and pipette lysate on membrane of Spin Filter. Close lid and centrifuge Receiver Tube with Spin Filter in a standard table top centrifuge for 1 min at 13,000 x g.
4. Open lid and add 300 μ l Wash Buffer I (see Preparation) to membrane of Spin Filter. Close lid and centrifuge for 30 sec at 13,000 x g.
5. Open lid and discard Receiver Tube with filtrate. Place Spin Filter in a new, labeled 2.0 ml Receiver Tube and add 750 μ l Wash Buffer II (see Preparation). Close lid and centrifuge for 30 sec at 13,000 x g.
6. Open lid and discard filtrate. Place Spin Filter in the same Receiver Tube, close lid and centrifuge for 2 min at 13,000 x g to eliminate residual ethanol.
7. Transfer Spin Filter to a new, labeled 1.5 ml Receiver tube (white). Add 200 μ l prewarmed Elution Buffer D to membrane of Spin Filter, close lid and incubate at room temperature for 1 min.
8. To elute DNA, centrifuge for 1 min at 6,000 x g.

DNA Isolation from Eukaryotic Cell Pellets

Best results with cell pellets are obtained from fresh cells. For long time storage, freeze material at -20°C to -80°C immediately after culture-harvesting. Avoid repeated freezing and thawing since this might accelerate degradation process. Inappropriate storage of sample material leads to decreased DNA quantity and quality. The DNA amount that can be purified with this kit from 10²-10⁵ cells depends on the type and number of cells.

Perform initial steps as indicated in section Preparation.

1. Pellet 10²-10⁵ eukaryotic cells in a labeled 1.5 ml reaction tube for 5 min at 300 x g. Remove supernatant carefully, but completely. Add 100 μ l Lysis Buffer M and 10 μ l Proteinase K (see Preparation) and vortex reaction tube thoroughly for at least 10 sec.
2. Incubate sample for 10 min at 56°C and 800 rpm in a thermal mixer.
3. Add 100 μ l Binding Buffer B6 and mix by pipetting up and down four times. Place a Spin Filter in a labeled 2.0 Receiver Tube (yellow) and pipette lysate on membrane of Spin Filter. Close lid and centrifuge Receiver Tube (yellow) with Spin Filter in a standard table top centrifuge for 1 min at 13,000 x g.
4. Open lid and add 300 μ l Wash Buffer I (see Preparation) to membrane of Spin Filter. Close lid and centrifuge for 30 sec at 13,000 x g.
5. Open lid and discard Receiver Tube with filtrate. Place Spin Filter in a new, labeled 2.0 ml Receiver Tube and add 750 μ l Wash Buffer II (see Preparation). Close lid and centrifuge for 30 sec at 13,000 x g.
6. Open lid and discard filtrate. Place Spin Filter in the same Receiver Tube, close lid and centrifuge for 2 min at 13,000 x g to eliminate residual ethanol.
7. Transfer Spin Filter to a new, labeled 1.5 ml Receiver tube (white). Add 200 μ l prewarmed Elution Buffer D to membrane of Spin Filter, close lid and incubate at room temperature for 1 min.
8. To elute DNA, centrifuge for 1 min at 6,000 x g.

DNA Isolation from Tissue

Best results with tissue (e.g. biopsy material, frozen section) are obtained from fresh specimens. For long time storage, freeze material at -20°C or -80°C immediately after withdrawal. Avoid repeated freezing and thawing since this might accelerate the degradation process. Inappropriate storage of sample material leads to decreased DNA quantity and quality. The DNA amount that can be purified with this kit from max. 5 mg tissue depends on the tissue type of the starting material. If applicable, thaw frozen section in Lysis Buffer M shortly before starting the isolation protocol.

Perform initial steps as indicated in section Preparation.

1. Transfer max. 5 mg tissue sample into a labeled 1.5 ml reaction tube. Grinding of the material, e.g. with a pipette tip, increases lysis efficiency. Add 100 μ l Lysis Buffer M and 10 μ l Proteinase K (see Preparation) and vortex reaction tube thoroughly for at least 10 sec.
2. Incubate sample for 30 min at 56°C and 800 rpm in a thermal mixer. If necessary, the lysis time can be increased.
3. Centrifuge sample for 2 min at maximum speed and transfer supernatant to a new, labeled 1.5 ml reaction tube (not provided).

4. Add 100 µl Binding Buffer B6 and mix by pipetting up and down four times. Place a Spin Filter in a labeled 2.0 Receiver Tube (yellow) and pipette lysate on membrane of Spin Filter. Close lid and centrifuge Receiver Tube with Spin Filter in a standard table top centrifuge for 1 min at 13,000 x g.
5. Open lid and add 300 µl Wash Buffer I (see Preparation) to membrane of Spin Filter. Close lid and centrifuge for 30 sec at 13,000 x g.
6. Open lid and discard Receiver Tube with filtrate. Place Spin Filter in a new, labeled 2.0 ml Receiver Tube and add 750 µl Wash Buffer II (see Preparation). Close lid and centrifuge for 30 sec at 13,000 x g.
7. Open lid and discard filtrate. Place Spin Filter in the same Receiver Tube, close lid and centrifuge for 2 min at 13,000 x g to eliminate residual ethanol.
8. Transfer Spin Filter to a new, labeled 1.5 ml Receiver tube (white). Add 200 µl prewarmed Elution Buffer D to membrane of Spin Filter, close lid and incubate at room temperature for 1 min.
9. To elute DNA, centrifuge for 1 min at 6,000 x g.

Limitations

The **GenoType® DNA Isolation Kit** is intended for use in *in vitro* diagnostic analysis, where DNA isolation is followed by an amplification or another enzymatic DNA modification.

Results generated with any diagnostic assay subsequent to an isolation with this test kit should be interpreted with regard to other clinical or laboratory findings.

Use of this assay is limited to qualified personnel well trained in the test procedure and familiar with molecular biological methods.

This kit must not be used for DNA isolation from stool samples. It was neither validated for isolation from bacteria, fungi, or viruses nor for isolation of RNA.

Appropriate sample storage (see chapter Test Procedure) is crucial for reproducible and high yields.

The performance evaluation of this kit was carried out using the protocols described above. The user is in charge of validating the applicability of the **GenoType® DNA Isolation Kit** for the intended specific use.

Troubleshooting

Low DNA yield due to

- a) insufficient cell lysis
 - Vortex sample thoroughly after addition of Lysis Buffer M and Proteinase K.
 - Extend incubation time with Lysis Buffer M and Proteinase K.
 - Cut tissue into small pieces. Extend lysis time until only few and small particles are left, (hair and bones, however, will not be lysed at all).
 - Ensure continuous shaking during lysis.
 - Avoid loss of activity of Proteinase K due to wrong storage or dilution (see chapter Test Procedure).
 - Use less starting material.
- b) inefficient binding of DNA to membrane
 - Use the correct amount of Binding Buffer B6.
 - Mix sample thoroughly prior to transfer to the Spin Filter.
- c) incomplete DNA elution
 - Extend centrifugation time for ethanol removal.
 - Extend incubation time with prewarmed Elution Buffer D to 5 min
 - Prewarm Elution Buffer D to 80°C.
 - Use larger volume of Elution Buffer D.
 - Elute twice each with 100 µl of Elution Buffer D.
- d) low DNA concentration
 - Elute DNA with lower volume of Elution Buffer D.
 - Ensure that Wash Buffer I and II were prepared correctly and used in the correct order (see chapter Test Procedure).
- e) too few cells in starting material
 - Ensure appropriate starting material.
 - Use more starting material (max. 50 µl whole blood or body fluids / max. 5 mg tissue / max. 10⁵ cells).
- f) incomplete flow-through of lysate through membrane
 - Extend centrifugation time until all the lysate has passed through the membrane.
 - Clogged membrane due to more than 10⁵ cells in starting material; reduce number of cells.
 - Clogged membrane due to cryoprecipitates; avoid repeated freezing and thawing, if cryoprecipitates appear, avoid to transfer them to the membrane of the Spin Filter.

Degradated DNA due to

- a) inappropriate storage of starting material
 - Freeze sample directly after withdrawal/harvest in liquid nitrogen and store at -80°C.
 - Avoid repeated freezing and thawing (aliquot sample, if necessary).
- b) old starting material
 - Provide new sample.

Problems in subsequent applications (e.g. PCR) due to

- a) ethanol in eluate
 - Ensure that centrifugation time and speed was followed as indicated in the protocol.
 - If necessary, extend the centrifugation time for ethanol removal.
- b) salt in eluate
 - Ensure that the Wash Buffers are equilibrated to room temperature.
 - Check Wash Buffers for salt precipitates. If necessary, warm carefully to up to max. 30°C to dissolve precipitates.
- c) use of too much/too little DNA
 - Optimize the amount of DNA used in the respective application.

Material Required but not Provided

- Adjustable pipettes for 20, 200, and 1000 µl
- Disposable gloves
- Disposable sterile pipette tips with filter
- ethanol p.a. 99,8%
- Graduated cylinder (250 ml)
- Reaction tubes (1.5 ml)
- Standard table top centrifuge
- Thermal mixer (precision +/-1°C)
- Vortexer
- Water (dd and molecular biology grade, if necessary)

Kit Contents

	supplied	
Spin Filters	12	48
Receiver Tubes 2.0 ml	24	2x 48
Receiver Tubes 1.5 ml	12	48
Lysis Buffer M	2 ml ready to use	3x 2 ml ready to use
Binding Buffer B6	2x 2 ml ready to use	15 ml ready to use
Proteinase K Lyophilisate	for 250 µl working solution	for 550 µl working solution
Wash Buffer I	30 ml (final volume 60 ml)	30 ml (final volume 60 ml)
Wash Buffer II	18 ml (final volume 60 ml)	18 ml (final volume 60 ml)
Elution Buffer D	2x 2 ml ready to use	15 ml ready to use
manual	1	1

Technical Data

Starting material: 1-50 µl whole blood (EDTA or citrate) / body fluids, 0.5-5 mg tissue, or 10²-10⁶ eukaryotic cells
 Yield: up to 2 µg genomic DNA, depending on nature and amount of starting material
 DNA purity: A₂₆₀/A₂₈₀ 1.7-2.0
 Assay time: 10-20 minutes (after lysis)

Hain Lifescience GmbH

Hardwiesenstraße 1, 72147 Nehren, Germany
<http://www.hain-lifescience.de>

51512

12

51548

48