Isolierung bakterieller DNA zum Nachweis parodontopathogener Markerkeime mittels High pure PCR template preparation Kit (Fa. Roche Diagnostics)

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12/2003

Isolation of DNA from Periodontopathogenic Bacteria with High Pure PCR Template Preparation Kit (Roche Diagnostics)

The working area must be free from amplified DNA. Before using for the first time, add the appropriate amounts of ethanol to Washing Buffer and Inhibitor Removal Buffer (see bottles). Heat two heating blocks to 72°C and 95°C, respectively. Prewarm 200-400 µl of Elution Buffer (see step 13 below) per sample to 72°C in a screw cap tube.

- 1) Per sample combine 180 µl of Tissue Lysis Buffer and 40 µl of Proteinase K solution in a tube and mix.
- 2) Add 220 µl of the above solution to each sample tube containing paper point/s. Vortex for 30 s.
- 3) Incubate for 10 min at 72°C. Briefly spin down.
- 4) Add 200 µl Binding Buffer, vortex for 15 s, and incubate for 5 min at 95°C. Briefly spin down.
- 5) Add 100 µl isopropanol, vortex for 15 s, and briefly spin down.
- 6) Insert labelled filter tube in collection tube and transfer solution from step 5 (without paper point/s) into the reservoir of the filter tube.
- 7) Close filter tube and spin for 1 min at 6000 x g (ap. 8000 rpm) in a standard table top centrifuge.
- 8) *optional*: Discard flow-through, re-insert filter tube in collection tube, add 500 µl Inhibition Removal Buffer and spin down (see step 7).
- 9) Discard flow-through, re-insert filter tube in collection tube, add 500 µl Washing Buffer and spin down (see step 7).
- 10) Repeat step 9.
- 11) *optional*: Insert filter tube in a new collection tube and spin down at maximum speed for 1 min.
- 12) Discard collection tube and insert filter tube in a clean and labelled 1.5 ml microcentrifuge tube.
- 13) For elution of DNA, add prewarmed Elution Buffer to filter tube. Use 200 µl when one or two paper points were present at the beginning and 400 µl when three or four points were present. Incubate for 1 min at room temperature and spin down for 1 min at 6000 x g. Discard spin column; the DNA is now in the eluate.
- 14) *optional*: In case a white precipitate is present in the eluate, spin samples again for 1 min at maximum speed and transfer supernatant to a clean and labelled tube.
- 15) Use 5 µl of the eluate for PCR.

Isolation of DNA from Periodontopathogenic Bacteria with QIAamp® DNA Mini Kit (QIAGEN)

The working area must be free from amplified DNA. Before using for the first time, add the appropriate amounts of ethanol to Buffers AW1 and AW2 (see bottles). If a precipitate has formed in Buffer ATL, dissolve in warm water. Heat two heating blocks to 70°C and 95°C, respectively.

- 1) Per sample combine 180 µl of Tissue Lysis Buffer and 20 µl of Proteinase K solution in a tube and mix.
- 2) Add 200 µl of the above solution to each sample tube containing paper point/s. Vortex for 30 s.
- 3) Incubate for 10 min at 70°C. Briefly spin down.
- 4) Add 200 µl Buffer AL, vortex for 15 s, and incubate for 5 min at 95°C. Briefly spin down.
- 5) Add 200 µl ethanol (96–100%), vortex for 15 s, and briefly spin down.
- 6) Insert labelled filter tube in collection tube and transfer solution from step 5 into the reservoir of the filter tube.
- 7) Close spin column and spin for 1 min at 6000 x g (ap. 8000 rpm) in a standard table top centrifuge.
- 8) Discard flow-through, re-insert spin column into collection tube, add 500 µl of Buffer AW1 and centrifuge as outlined in step 7.
- 9) Discard flow-through, re-insert spin column into collection tube, add 500 µl of Buffer AW2 and centrifuge for 3 min at maximum speed.
- 10) *optional*: Discard collection tube, insert spin column in a clean collection tube and spin for 1 min at maximum speed.
- 11) Discard collection tube and insert filter tube in a clean and labelled 1.5 ml microcentrifuge tube.
- 12) Open spin column and add Buffer AE. Use 200 µl when one or two paper points were present at the beginning and 400 µl when three or four points were present. Incubate for 1 min at room temperature and spin for 1 min at 6000 x g. Discard spin column; the DNA is now in the eluate.
- 13) *optional*: In case a white precipitate is present in the eluate, spin samples again for 1 min at maximum speed and transfer supernatant to a clean and labelled tube.
- 14) Use 5 µl of the eluate for PCR.