Starterset Stanze

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Starter Set Punch

Punch and Punching Mats for GENO•CARD®

The punch and the punching mats allow to cut discs from samples previously spotted on a **GENO•CARD®**. Subsequently, the discs are used for amplification in test kits from Hain Lifescience validated for this DNA isolation technique. Compatible kits have an according note in the manual.

The punch has a sharp cutting tip for punching discs of 1.2 mm in diameter and a plunger to eject the discs. The punching mat is made from rigid plastic and is placed inside the **GENO•CARD®** for punching.

Storage and Precautions

Store punch and punching mat protected from light at a dry place at room temperature. Due to the sharp cutting tip and the risk of injury the punch should be handled with care. Put on protective cap after use.

To avoid contaminations, the following cleaning procedures should be applied at regular intervals:

To clean the punch, screw off the tip, rinse it with 1.5% sodium hypochlorite solution and subsequently with demineralized water. Perform the same procedure with plunger. Dry with compressed air or by rinsing with 99.8% ethanol. Clean outer surface of punch with 70% ethanol.

After use, remove punching mat from **GENO•CARD®**, clean thoroughly with 1.5% sodium hypochlorite solution, rinse with demineralized water and air dry.

Collecting a Disc

A detailed description of the collection of discs from sample fields can be found in the manual for the **GENO•CARD®**.

If clear signs of usage are visible on the punching mat, it should be replaced. Please note that each mat features four areas for punching. If the punch tip has become blunt, a spare tip can be ordered from Hain Lifescience (punch tip incl. 5 punching mats; order no. PR027).

Limitations

The performance evaluation of the **GENO•CARD®** DNA collection system was carried out with compatible test kits from Hain Lifescience only, applying the conditions indicated in the respective manual. The user is in charge of validating the applicability of the **GENO•CARD®** for any intended specific use.

Contents

punch	1
punching mats	5
manual	1

Ordering Information

GENO•CARD® for 3x 4 samples	order no. G001
GENO•CARD® for 24x 4 samples	order no. G00196
Spare Tip Set (punch tip incl. 5 punching mats)	order no. PR027

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Hain Lifescience GmbH