GenoType[®] MRSA

VER 2.0

Deutsch: S. 2-12

English: p. 13-22

Français : p. 23-33

Italiano: p. 34-44

Español: p. 45-55

Português: p. 56-66

Česky: s. 67-77

08/2008

HAIN
LIFESCIENCE

GenoType® MRSA

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GenoType® MRSA Molecular Genetic Assay for Fast Identification of Methicillin-Resistant Staphylococci from Cultured Material

Methodology

The **GenoType® MRSA** test is based on the **DNA•STRIP®** technology and permits the identification of *S. aureus* and *S. epidermidis* strains. At the same time the methicillin resistance mediating *mecA* gene and the bicomponent cytotoxin virulence factor PVL (Panton-Valentine leukocidin) can be detected. The whole procedure is divided into three steps: DNA extraction from cultured material (bacteria freshly grown on culture plates; the necessary reagents are not provided), a multiplex amplification with biotinylated primers (the necessary thermostable DNA polymerase is not provided), and a reverse hybridization. The hybridization includes the following steps: chemical denaturation of the amplification products, hybridization of the single-stranded, biotin-labeled amplicons to membrane-bound probes, stringent washing, addition of a streptavidin/alkaline phosphatase (AP) conjugate, and an AP mediated staining reaction. A template ensures the easy and fast interpretation of the banding pattern obtained.

Storage and Precautions

Store Primer/Nucleotide Mix (PNM) at 2-8°C upon arrival isolated from any potential source of contaminating DNA. If longer storage (more than 4 weeks) is required, store at -20°C. In order to avoid repeated freezing and thawing, aliquot PNM. Store all other kit components at 2-8°C. Do not use the reagents beyond their expiry date. Patient specimens and cultures made from patient specimens must always be considered as potentially infectious. Samples from risk patients and cultures made from those samples must always be labeled and handled under suitable safety conditions. Observe all federal, state, and local safety and environmental regulations. Always wear suitable protective clothing and gloves.

Observe the usual precautions for amplification set-up. It is essential that all reagents and materials used for DNA extraction and amplification set-up are free from DNases.

When handling kit reagents, the following special safety measures must be applied:

The **Denaturation Solution** (DEN) contains <2% NaOH and is irritating to eyes and skin (R 36/38 and S 26-37/39-45).

The **Substrate Concentrate** (SUB-C) contains Dimethyl Sulfoxide and is irritating (R 36/37/38, S 23-26-36).

For additional information, please refer to material safety data sheets which can be downloaded from: www.hain-lifescience.com/products/msds.html

Quality Control

In order to validate the correct performance of the test and the proper functioning of reagents, each strip includes 2 control zones:

- a Conjugate Control zone to check the binding of the conjugate on the strip and a correct chromogenic reaction
- a Universal Control zone which detects, as known, all bacterial species

DNA Extraction

The preferred starting material for DNA extraction consists of bacteria freshly grown on culture plates. The working area must be free from amplified DNA. When using bacteria grown in liquid medium, the probability of examining mixed populations is higher thus complicating interpretation of results. Any DNA extraction procedure producing amplifiable DNA from bacteria can be used. The following quick protocol normally also yields DNA suitable for amplification:

- 1a. When using bacteria grown on solid media, collect up to 5 colonies with an inoculation loop and suspend in 150 µl of water.
- 1b. When using bacteria grown in liquid media, directly apply 1 ml. Pellet bacteria by spinning for 5 min in a standard table top centrifuge at 5000-6000 x g (ap. 7500-8000 rpm). Discard supernatant, resuspend bacteria in 1 ml of water by vortexing and spin down as before. Discard supernatant and resuspend bacteria in 50 µl of water.
Using too much bacterial material for extraction may lead to inhibition of PCR. When densely grown cultures are used, the amount of starting material should be reduced.
2. Incubate solution from 1a or 1b at 95°C for 15 min.
3. Incubate for 15 min in an ultrasonic bath.

4. Spin down for 5 min at maximum speed and directly use 5 μl of the supernatant for PCR. In case DNA solution is to be stored for an extended time period, transfer supernatant to a new tube.

Amplification

Prepare the amplification mix (45 μl) in a DNA-free room. The DNA sample should be added in a separate area.

Per tube mix:

- 35 μl PNM
- 5 μl 10x polymerase incubation buffer – not provided
- x μl MgCl_2 solution¹⁾ – not provided
- 1-2 unit(s) thermostable DNA polymerase (refer to manual) – not provided
- y μl water to obtain a volume of 45 μl (not considering volume of enzyme) – not provided
- Add 5 μl freshly extracted DNA solution (20-100 ng DNA) leading to a final volume of 50 μl (not considering volume of enzyme).

¹⁾ Depending on the enzyme/buffer system used, the optimal MgCl_2 concentration may vary between 1.5 and 2.5 mM. Please note that some incubation buffers already contain MgCl_2 .

Determine the number of samples to be amplified (number of samples to be analyzed plus control samples). A negative control sample, for example, contains 5 μl of water instead of DNA solution. Prepare a master mix containing all reagents except for DNA solution and mix well (do not vortex). Aliquot 45 μl in each of the prepared PCR tubes.

Amplification profile:

15 min²⁾ 95°C 1 cycle

20 sec 95°C
30 sec 60°C } 22 cycles

²⁾ When using certain hot start DNA polymerases, the time interval of this step has to be reduced (please refer to manual of the enzyme).

Amplification products can be stored at +4 to -20°C.

For checking the amplification reaction, 5 µl of each sample might be directly applied to a gel without the addition of loading buffer. The amplicons have a length of 63 bp (Universal Control fragment), 71 bp (*mecA* fragment), 77 bp (PVL fragment), 85 bp (*S. aureus*-specific fragment), and 125 bp (*S. epidermidis*-specific fragment), respectively.

Hybridization

Preparation

Prewarm shaking water bath/**TwinCubator®** to **45°C**; the maximum tolerated deviation from the target temperature is +/-1°C. Prewarm solutions HYB and STR to 37-45°C before use. The reagents must be free from precipitates (note, however, that solution CON-D is opaque). Mix if necessary. Warm the remaining reagents with the exception of CON-C and SUB-C to room temperature. Using a suitable tube, dilute Conjugate Concentrate (CON-C, orange) and Substrate Concentrate (SUB-C, yellow) 1:100 with the respective buffer (**CON-C with CON-D, SUB-C with SUB-D**) in the amounts needed. Mix well and bring to room temperature. For each strip, add 10 µl concentrate to 1 ml of the respective buffer. Dilute CON-C before each use. Diluted SUB-C is stable for 4 weeks if stored at room temperature and protected from light.

Note: This test can also be performed using a shorter protocol which can be downloaded from: www.hain-lifescience.de/pdf/mrsa_short.pdf

1. **Dispense 20 µl of Denaturation Solution (DEN, blue) in a corner of each of the wells used.**
2. **Add to the solution 5-20 µl of amplified sample, pipette up and down to mix well and incubate at room temperature for 5 minutes.**

Meanwhile, take strips out of the tube using tweezers and mark them with a pencil underneath the colored marker. Always wear gloves when handling strips.

3. **Carefully add to each well 1 ml of prewarmed Hybridization Buffer (HYB, green). Gently shake the tray until the solution has a homogenous color.**

Take care not to spill solution into the neighboring wells.

4. Place a strip in each well.

The strips must be completely covered by the solution and the coated side (identifiable by the colored marker near the lower end) must face upward. Using tweezers, turn over strips which might have turned when immersed in the solution. Carefully clean tweezers after each use to avoid contamination. This also applies to all following steps.

5. Place tray in shaking water bath/TwinCubator® and incubate for 30 minutes at 45°C.

Adjust the shaking frequency of the water bath to achieve a constant and thorough mixing of the solution. To allow adequate heat transfer, the tray must be dipped into the water to at least 1/3 of its height.

6. Completely aspirate Hybridization Buffer.

For example, use a Pasteur pipette connected to a vacuum pump.

7. Add 1 ml of Stringent Wash Solution (STR, red) to each strip and incubate for 15 minutes at 45°C in shaking water bath/TwinCubator®.

8. Work at room temperature from this step forward.

Completely remove Stringent Wash Solution.

Pour out Wash Solution in a waste container and remove all remaining fluid by turning the tray upside down and gently striking it on an absorbent paper. This also applies to all other wash steps.

9. Wash each strip once with 1 ml of Rinse Solution (RIN) for 1 minute on shaking platform/TwinCubator® (pour out RIN after incubation).

10. Add 1 ml of diluted Conjugate (see above) to each strip and incubate for 30 minutes on shaking platform/TwinCubator®.

11. Remove solution and wash each strip twice for 1 minute with 1 ml of Rinse Solution (RIN) and once for 1 minute with approx. 1 ml of distilled water (e. g. use wash bottle) on shaking platform/TwinCubator® (pour out solution each time). Make sure to remove any trace of water after the last wash.

12. Add 1 ml of diluted substrate (see above) to each strip and incubate protected from light without shaking. Depending on the test conditions (e. g. room temperature), the substrate incubation time can vary between 3 and 20 minutes. Extended substrate incubation times can lead to increased background staining and might impair interpretation of the results.

13. Stop reaction by briefly rinsing twice with distilled water.


14. Using tweezers, remove strips from the tray and dry them between two layers of absorbent paper.

Evaluation and Interpretation of Results

Paste strips and store protected from light. An evaluation sheet is provided with the kit and can be downloaded from:

www.hain-lifescience.de/pdf/mrsa_evaluation.pdf

The supplied template serves as an aid for evaluation and is aligned with the Conjugate Control band of the strip. Each strip has a total of 6 reaction zones (see figure).

Note: The strip is not displayed in original size and must not be used for interpretation purposes.

Conjugate Control (CC)

A line must develop in this zone, documenting the efficiency of conjugate binding and substrate reaction.

Universal Control (UC)

This probe is directed against a highly conserved DNA region appearing in all bacterial species tested so far. It indicates the presence of bacterial DNA and the correct implementation of DNA extraction and amplification. Since the source material for this test is a bacterial culture, this reaction zone should be developed. The amplification of the

control sequence can, however, be suppressed by the amplification of the specific sequence(s).

Please be aware that the Universal Control is not a PCR control, i.e. when implementing a contamination control (water instead of DNA solution in the PCR preparation) this reaction zone remains negative.

Other bands

Specific probes; for evaluation see figure/template.

In case none of the other bands show a positive signal, the control zones have to stain positive. If this is not the case a false negative reaction occurred and the test has to be repeated.

Not all bands of a strip have to show the same signal strength.

Limitations

Prior to amplification, DNA has to be extracted from cultured bacteria using a suitable method. It must be ensured that the template DNA is efficiently amplified during the amplification reaction.

The test only works within the limits of the genomic regions the primers and probes were chosen from. Potential sequence analysis remains reserved to resuming investigations.

The presence of multiple bacterial species in the sample to be analyzed might hamper the interpretation of the test.

As with any detection system on hybridization basis the test system on hand bears the possibility that sequence variations in the genomic regions the primers and probes were chosen from but the detection of which the test was not designed for may lead to false results. Due to the high variability of bacterial genomes it is possible that certain sub-types might not be detected. The test reflects the state of knowledge of Hain Lifescience.

Use of this assay is limited to qualified personnel well trained in the test procedure and familiar with molecular biological methods.

Performance evaluation of this assay was carried out using the HotStarTaq polymerase from Qiagen (Hilden, Germany).

Troubleshooting

Overall weak or no signals (including Conjugate Control zone)

- Room temperature too low or reagents not equilibrated to room temperature.
- No or too little amount of CON-C and/or SUB-C used.

Weak or no signals except for Conjugate Control zone

- Quality and/or quantity of extracted DNA do not allow an efficient amplification. Check amplicon on a gel. In case no amplicon is visible, repeat DNA extraction and amplification. If necessary, try a different DNA extraction method (see chapter DNA Extraction). Please note that the amplification can also be inhibited by using too much bacterial material.
- Incubation temperature too high.

No homogeneous staining

- Strips were not completely immersed during incubation steps.
- Tray was not shaken properly.

High background color

- CON-C and/or SUB-C used too concentrated.
- Washing steps were not performed with the necessary care.
- Wash solutions too cold.

Unexpected result

- Wrong incubation temperature.
- Hybridization Buffer and/or Stringent Wash Solution were not properly prewarmed or mixed.
- Contamination of extracted DNA and/or amplification agents with extracted and/or amplified DNA. In case amplification agents are contaminated a negative control sample also shows the respective banding pattern.
- Contamination of neighboring wells by spillage during addition of Hybridization Buffer.
- Depending on the amount of amplified DNA used and on the specific reaction conditions, a strong and fast color development may occur. In such cases, dis-

- continue the substrate incubation as soon as the signals are clearly visible in order to prevent the development of cross-hybridizing bands.
- No pure culture as starting material.
 - The bacterial species present in the sample can not be detected with this test.

Material Required but not Provided

- Absorbent paper
- Adjustable pipettes for 10, 20, 200, and 1000 μ l
- Calibrated thermometer
- DNA extraction reagents for amplification use as well as necessary equipment
- Disposable gloves
- Disposable sterile pipette tips with filter
- Graduated cylinder
- PCR tubes, DNase and RNase free
- Shaking water bath/**TwinCubator**[®]
- Shaking platform/**TwinCubator**[®]
- Thermal cycler (heating rate: 3°C/sec, cooling rate: 2°C/sec, precision: +/-0.2°C)
- Thermostable DNA polymerase with buffer (recommendation: hot start enzyme, extension rate: 2-4 kb/min at 72°C, half-life: 10 min at 97°C, 60 min at 94°C, amplification efficiency: >10⁵ fold)
- Timer
- Tweezers
- Water (molecular biology grade)

Kit Contents

	Supplied	
Membrane strips coated with specific probes (STRIPS)	12	96
Primer Nucleotide Mix (PNM) contains specific primers, nucleotides, dye	0.5 ml	4 ml
Denaturation Solution (DEN) ready to use contains <2% NaOH, dye	0.3 ml	2.4 ml
Hybridization Buffer (HYB) ready to use contains 8-10% anionic tenside, dye	20 ml	120 ml
Stringent Wash Solution (STR) ready to use contains >25% of a quaternary ammonium compound, <1% anionic tenside, dye	20 ml	120 ml
Rinse Solution (RIN) ready to use contains buffer, <1% NaCl, <1% anionic tenside	50 ml	360 ml
Conjugate Concentrate (CON-C) concentrate contains streptavidin-conjugated alkaline phosphatase, dye	0.2 ml	1.2 ml
Conjugate Buffer (CON-D) contains buffer, 1% blocking reagent, <1% NaCl	20 ml	120 ml
Substrate Concentrate (SUB-C) concentrate contains Dimethyl Sulfoxide, substrate solution	0.2 ml	1.2 ml
Substrate Buffer (SUB-D) contains buffer, <1% MgCl ₂ , <1% NaCl	20 ml	120 ml
tray, evaluation sheet	1 of each	4 of each
manual, template	1 of each	1 of each

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Hain Lifescience GmbH