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07/2007

DNA Isolation from Mycobacteria using the QIAamp DNA Mini Kit from Qiagen

BE AWARE OF RISK OF INFECTION WHEN HANDLING MATERIAL CONTAINING MYCOBACTERIA!!!

The working area must be free from amplified DNA. Before using for the first time, add the appropriate amounts of ethanol to Buffers AW1 and AW2 (see bottles). If a precipitate has formed in Buffer ATL, dissolve in warm water. Heat two heating blocks to 70°C and 95°C, respectively. It is crucial to heat samples to 95°C for at least 20 min in order to inactivate vegetative bacteria.

- 1) **a** When using bacteria grown on solid media, collect an inoculation loop of bacteria and suspend in 300 µl water (molecular biology grade).
- b** When using bacteria grown in liquid media, directly apply 1 ml.

Using too many bacteria for isolation may lead to inhibition of PCR. When densely grown cultures are used, the amount of starting material should be reduced.

- 2) Incubate bacteria from 1a or 1b at 95°C for 20 min.
- 3) Pellet bacteria by spinning for 15 min in a standard table top centrifuge at 10 000 x g (ca 10 000-11 000 rpm).
- 4) Discard supernatant, resuspend bacterial pellet in 180 µl Buffer ATL, add 20 µl Proteinase K, vortex briefly, and incubate at 70°C for 10 min.
- 5) Briefly spin down, add 200 µl Buffer AL, and vortex for 15 s.
- 6) Briefly spin down, add 200 µl ethanol (96-100%), and vortex for 15 s.
- 7) Insert labelled filter tube in collection tube and transfer solution from step 6 into the reservoir of the filter tube.
- 8) Close filter tube and spin for 1 min at 6000 x g (ca 8000 rpm).
- 9) Discard flow-through, re-insert spin column into collection tube, add 500 µl of Buffer AW1 and centrifuge as outlined in step 8.
- 10) Discard flow-through, re-insert spin column into collection tube, add 500 µl of Buffer AW2 and centrifuge for 3 min at maximum speed.
- 11) optional: Discard collection tube, insert spin column in a clean collection tube and spin for 1 min at maximum speed.
- 12) Discard collection tube and insert filter tube in a clean and labelled 1.5 ml microcentrifuge tube.
- 13) For elution of DNA, add 100 µl Buffer AE to filter tube, incubate for 1 min at room temperature and spin down (1 min, 6000 x g). Discard spin column; the DNA is now in the eluate.

- 14) optional: In case a white precipitate is present in the eluate, spin samples again for 1 min at maximum speed and transfer supernatant to a clean and labelled tube.
- 15) Use 5 μ l of the eluate for PCR.

07/2007

DNA Isolation from Mycobacteria using the High Pure PCR Template Preparation Kit from Roche Diagnostics

BE AWARE OF RISK OF INFECTION WHEN HANDLING MATERIAL CONTAINING MYCOBACTERIA!!!

The working area must be free from amplified DNA. Before using for the first time, add the appropriate amounts of ethanol to Washing Buffer and Inhibitor Removal Buffer (see bottles). Heat two heating blocks to 72°C and 95°C, respectively. Prewarm ca 120 µl of Elution Buffer per sample to 72°C in a screw cap tube. It is crucial to heat samples to 95°C for at least 20 min in order to inactivate vegetative bacteria.

- 1) **a** When using bacteria grown on solid media, collect an inoculation loop of bacteria and suspend in 300 µl water (molecular biology grade).
- b** When using bacteria grown in liquid media, directly apply 1 ml.

Using too many bacteria for isolation may lead to inhibition of PCR. When densely grown cultures are used, the amount of starting material should be reduced.

- 2) Incubate bacteria from 1a or 1b at 95°C for 20 min.
- 3) Pellet bacteria by spinning for 15 min in a standard table top centrifuge at 10 000 x g (ca 10 000-11 000 rpm).
- 4) Discard supernatant, resuspend bacterial pellet in 200 µl Tissue Lysis Buffer, add 40 µl Proteinase K, briefly vortex, and incubate at 72°C for 10 min.
- 5) Briefly spin down, add 200 µl Binding Buffer, and vortex for 15 s.
- 6) Briefly spin down, add 100 µl isopropanol, and vortex for 15 s.
- 7) Insert labelled filter tube in collection tube and transfer solution from step 6 into the reservoir of the filter tube.
- 8) Close filter tube and spin for 1 min at 6000 x g (ca 8000 rpm).
- 9) *optional*: Discard flow-through, re-insert filter tube in collection tube, add 500 µl Inhibition Removal Buffer and spin down (see step 8).
- 10) Discard flow-through, re-insert filter tube in collection tube, add 500 µl Washing Buffer and spin down (see step 8).
- 11) Repeat step 10.
- 12) *optional*: Insert filter tube in a new collection tube and spin down at maximum speed for 1 min.
- 13) Discard collection tube and insert filter tube in a clean and labelled 1.5 ml microcentrifuge tube.

- 14) For elution of DNA, add 100 μ l prewarmed Elution Buffer to filter tube, incubate for 1 min at room temperature and spin down (1 min, 6000 x g). Discard spin column; the DNA is now in the eluate.
- 15) *optional*: In case a white precipitate is present in the eluate, spin sample again for 1 min at maximum speed and transfer supernatant to a clean and labelled tube.
- 16) Use 5 μ l of the eluate for PCR.

07/2007

Quick Method for DNA Isolation from Mycobacteria

BE AWARE OF RISK OF INFECTION WHEN HANDLING MATERIAL CONTAINING MYCOBACTERIA!!!

The working area must be free from amplified DNA. Heat a heating block or water bath to 95°C. Provide an ultrasonic bath. It is crucial to heat samples to 95°C for at least 20 min in order to inactivate vegetative bacteria.

Liquid medium

- 1) When using bacteria grown in liquid media, directly apply 1 ml.
Using too many bacteria for isolation may lead to inhibition of PCR. When densely grown cultures are used, the amount of starting material should be reduced.
- 2) Pellet bacteria from 1 by spinning for 15 min in a standard table top centrifuge at 10 000 x g (ca 10 000-11 000 rpm).
- 3) Discard supernatant, resuspend bacteria in 100-300 µl of distilled water (see above) and incubate for 20 min at 95°C (heating block or boiling water bath).
- 4) Incubate for 15 min in an ultrasonic bath.
- 5) Spin down for 5 min at maximum speed and directly use 5 µl of the supernatant for PCR. In case DNA solution is to be stored for an extended time period, transfer supernatant to a clean and labelled tube.

Solid medium

- 1) When using bacteria grown on solid medium, collect an inoculation loop of bacteria and suspend in 300 µl water (molecular biology grade).
Using too many bacteria for isolation may lead to inhibition of PCR. When densely grown cultures are used, the amount of starting material should be reduced.
- 2) Incubate bacteria suspension for 20 min at 95°C (heating block or boiling water bath).
- 3) Incubate for 15 min in an ultrasonic bath.
- 4) Spin down for 5 min at maximum speed and directly use 5 µl of the supernatant for PCR. In case DNA solution is to be stored for an extended time period, transfer supernatant to a clean and labelled tube.

07/2007.